


Assay for (1→3)-β-D-Glucan in Serum

FUNGITELL STAT®

Instructions For Use



ASSOCIATES OF
CAPE COD
INCORPORATED

124 Bernard E. Saint Jean Drive • E. Falmouth, MA 02536 USA

Telephone: (508) 540-3444
Toll-Free: (888) 395-2221
Fax: (508) 540-8680
Technical Support: (800) 848-3248
Customer Service: (800) 525-8378





PN002603-en Rev5 2023-06-13

Visit www.acciusa.com for instructions for use in your language.

This product is for In Vitro Diagnostic Use and Professional Use only.

1. Intended Use
The Fungitell STAT® assay is a protease zymogen-based colorimetric assay for the qualitative detection of (1→3)-β-D-glucan in the serum of patients with symptoms of, or medical conditions predisposing the patient to, invasive fungal infection. The serum concentration of (1→3)-β-D-glucan, a major cell-wall component of various medically important fungi¹, can be used as an aid in the diagnosis of deep-seated mycoses and fungemias². A positive result does not indicate which genus of fungi may be causing infection.

(1→3)-β-D-glucan index values should be used in conjunction with other diagnostic procedures, such as microbiological culture, histological examination of biopsy samples and radiological examination.

2. Summary and Explanation
There is an increasing incidence of fungal infections by opportunistic pathogens, especially in immunocompromised patients^{3,4,5}. Invasive fungal diseases, as opportunistic infections, are common among hematological malignancy and AIDS patients and account for a growing number of nosocomial infections, particularly among organ transplant recipients and other patients receiving immunosuppressive treatments^{6,7}. Many fungal diseases are acquired by inhaling fungal spores originating from the soil, plant detritus, air-handling systems and/or exposed surfaces. Some opportunistic fungi are present in/on human skin, the intestinal tract, and mucous membranes^{8,9}. Diagnosis of invasive mycoses and fungemias is usually based on non-specific diagnostic or radiological techniques. Recently, biological markers of fungal infection have been added to the available diagnostic methods².

Opportunistic fungal pathogens include *Candida spp.*, *Aspergillus spp.*, *Fusarium spp.*, *Trichosporon spp.*, *Saccharomyces cerevisiae*, *Acremonium spp.*, *Coccidioides immitis*, *Histoplasma capsulatum*, *Sporothrix schenckii*, *Exserohilum rostratum*, and *Pneumocystis jirovecii*. The (1→3)-β-D-glucan produced by these organisms, and others, can be detected by the Fungitell STAT® assay^{1,5,10,11}.

3. Principle of the Procedure
The Fungitell STAT® (cat# FT007, Associates of Cape Cod, Inc.) assay is a design modification to the Fungitell® (cat# FT001, Associates of Cape Cod, Inc. or ACC) assay format. The Fungitell STAT® assay (2019 CE-marked device) was developed to answer the need for a single use test format and smaller kit size relative to the 96-well plate format of the Fungitell® (USA predicate and 2008 CE-marked device) assay.

The Fungitell STAT® assay provides a qualitative measurement of (1→3)-β-D-glucan. The assay is based upon a modification of the *Limulus* Amebocyte Lysate (LAL) pathway^{12,13,14,15}. **Figure 1.** The Fungitell STAT® Reagent is modified to eliminate bacterial endotoxin reactivity and, thus, to only react to (1→3)-β-D-glucan, through the Factor G-mediated side of the pathway. (1→3)-β-D-glucan activates Factor G, a serine protease zymogen. The activated Factor G converts the inactive pro-clotting enzyme to the active clotting enzyme, which in turn cleaves the para-nitroanilide Boc-Leu-Gly-Arg-pNA, creating a chromophore, para-nitroaniline (pNA), that absorbs at 405 nm. The Fungitell STAT® kinetic assay, described below, is based upon the determination of the rate of optical density increase produced by a patient serum sample. This rate is compared to the rate of optical density increase of the Fungitell STAT® Standard to produce an index. The Fungitell STAT® Standard is calibrated at 80 +/- 8 pg/mL which is the Positive cut-off for the Fungitell® assay. This patient serum sample index value is qualitatively interpreted as a Negative, Indeterminate, or Positive result according to the index value ranges provided in **Table 1** below.




Figure 1. *Limulus* Amebocyte Lysate Pathway

Table 1. Fungitell STAT® Index Ranges	
Result	Index Value
Negative	≤ 0.74
Indeterminate	0.75 – 1.1
Positive	≥ 1.2

4. Materials Supplied with the Fungitell STAT® product
The Fungitell STAT® product is for *in vitro* diagnostic use.


The following materials supplied with each product are sufficient for a total of 10 reactions (based on the 10 tubes of Fungitell STAT® Reagent). Each product also contains 5 Fungitell STAT® Standard tubes.

- Fungitell STAT® Reagent, a lyophilized (1→3)-β-D-glucan specific LAL (10 tubes)
The Fungitell STAT® Reagent is composed of Limulus (i.e., horseshoe crab) Amebocyte Lysate, Boc-Leu-Gly-Arg-pNA colorimetric substrate and Tris buffer. It does not contain human or mammalian proteins. Fungitell STAT® Reagent is free of interfering levels of (1→3)-β-D-glucan.
- Fungitell STAT® Glucan Standard (5 tubes) lyophilized (1→3)-β-D-glucan.
*The Fungitell STAT® Glucan Standard is composed of D-lactose and (1→3)-β-D-glucan derived from *Saccharomyces cerevisiae* yeast extract.*
Internal control: The Fungitell STAT® Standard (1→3) β-D-glucan concentration is calibrated to the positive limit value of the Fungitell® product (USA predicate and CE mark 2008) and against an internal reference standard. The Fungitell STAT® Standard, contains a known amount of glucan. The resulting values are described in the Quality Control section and serve as internal control for the Fungitell STAT® assay.
- Instructions for Use
- Quick Visual Guide

5. Materials Required but not Supplied
All materials must be free of interfering glucan.

- LAL Reagent Water* (5.5 mL vial, catalog # W0051-10)
- Alkaline Pretreatment Solution 0.125 M KOH and 0.6 M KCl* (2.5 mL vial, catalog #APS51-5)
- Pipettes capable of delivering 20-200 µL and 100-1000 µL volumes
- Pipette tips* (250 µL catalog # PPT25 and 1000 µL catalog # PPT10)
- Long Pipette tips* (20-200 µL, catalog # TPT50)
- Test tubes* for patient sample preparation and combining serum pretreatment solution. (12 x 75 mm, catalog # TB240-5)
- Tube reader and kinetic assay software
 - PKF08 Incubating 8-Well Tube Reader (PKF08-1, Lab Kinetics, LLC)** with Beta Glucan Analytics (BG Analytics® or BG Analytics® Software), BG Analytics® Software Manual and BG Analytics® System Verification Protocol** (BGA007, Associates of Cape Cod, Inc.).
The PKF08 device and BG Analytics® software are supplied by Associates of Cape Cod, Inc. (catalog# PKF08-PKG**). The PKF08-PKG has been validated for use with the Fungitell STAT® test. Or...
 - Incubating (37°C) tube reader capable of reading at 405 nm and 495 nm with a range of at least 0 – 1.0 Absorbance Units, coupled with appropriate computer-based kinetic assay software capable of analyzing reaction kinetics as well as supporting the review of the criteria listed in the Quality Control section of the IFU.
- Sterile, glucan-free, tubes for aliquoting samples. Tubes that are certified to be RNase, DNase, and pyrogen-free can be used.
- Parafilm®

** These products, supplied by Associates of Cape Cod, Inc. (ACC), are certified free of interfering glucans. **User Manuals can be downloaded from ACC website: www.acciusa.com.*

- Reagent Storage**
 - Store the kit, as supplied, at 2-8°C in the dark.
 - Fungitell STAT® Reagent and Fungitell STAT® Standard are designed to be used up to 1 hour after reconstitution.
-  **Warnings and Precautions**
 - Do not pipette any material by mouth. Do not smoke, eat or drink in areas where specimens or kit reagents are handled.
 - Follow operational and local safety regulations.
 - Wear protective gloves when handling biological samples that may be infectious or dangerous. The gloved hands should be considered contaminated at all times; keep your gloved hands away from your eyes, mouth and nose. Wear an eye guard and surgical mask if there is a possibility of aerosol contamination.
 - Products with damaged contents should not be used.
 - Disposal: Residues of chemicals and preparations are generally considered to be hazardous wastes. The disposal of this type of waste is regulated by national and regional laws and regulations. Contact your local authorities or waste management companies for advice on the disposal of hazardous waste.
 - The Safety data sheets** for the Fungitell STAT® Reagent, Fungitell STAT® Standard, LAL Reagent Water and Alkaline Pretreatment Solution can be downloaded from ACC website: www.acciusa.com.

7.1 Procedural Precautions
The Fungitell STAT® assay requires rigorous attention to technique and the testing environment. Thorough training of the technician in the assay method and in the avoidance of contamination is critical for the effectiveness of the assay.

- Establish a clean environment in which to perform the assay.
- Note that glucan as well as fungal particles contamination from the human body, clothes, containers, water and airborne dust may cause interference with the Fungitell STAT® test.
- Possible sources of contamination include cellulose-containing materials such as gauze, paper wipes and cardboard, glass pipettes with cotton plugs and pipette tips with cellulose filters. Surgical gauze bindings and sponges can also secrete high amounts of (1→3)-β-D-glucan^{21,22}. For other patient-related sources of contamination, see the Limitations section of the test.
- Use the open vials with alkaline pretreatment solution and LAL reagent water immediately and if potential contamination is a concern, do not re-use these materials.
- The Fungitell STAT® Reagent and the Fungitell STAT® Standard are released as a paired batch. For this reason, no Fungitell STAT® Reagent and Fungitell STAT® Standard components from other product batches should be used. Therefore, it is recommended to dispose of any remaining Fungitell STAT® Standards as soon as all Fungitell STAT® Reagent tubes contained in a package have been used up.

- Do not use materials beyond their expiry date.

- Specimen Handling**
 - Blood collection and preparation of serum shall be carried out in accordance with applicable local regulations. Specimen Collection: Blood samples may be collected in sterile serum preparation tubes or serum separator tubes (SST) for the preparation of serum.
 - Specimen Storage: Serum samples can be stored at 2-8°C for up to 15 days, or frozen at -20°C for up to 27 days or -80°C for up to 4 years.
 - Specimen Labeling: Specimens should be clearly labeled according to the approved practices of the medical institution (laboratory).

- Notes on Testing:**
 - Use good laboratory practices according to your local regulations. This assay is sensitive to contamination and pipetting inaccuracy.
 - In order to ensure the safety of the operator while working with serum samples and to reduce the potential for contamination by (1→3)-β-D-glucan from the environment during the process, it is recommended to work in a biological safety cabinet.
 - To reduce unnecessary glass vial movements in and out of the biological safety cabinet, it is recommended to bring the vortex device within the biological safety cabinet (as long as the critical airflow is maintained).
 - It is recommended to use long pipette tips to help prevent cross-contamination between vials.
 - A Fungitell STAT® Standard (red cap and red line label) should always be processed under the same conditions and at the same time as the patient sample(s) within a run. This is critical since the outcome of the assay is an Index (sample/standard) of the kinetic reaction rates (or slopes, OD/sec) from the Patient sample and the Fungitell STAT® Standard.
 - It is recommended to use separate tube racks during the procedure, one for the sample preparation tubes and one for the reagent tubes. To avoid confusion and cross contamination.
 - It is recommended to place the Fungitell STAT® Standard at a defined and consistent position within the tube rack, incubator and reader. In the PKF08 Reader, use the first well on the left which is labeled “Standard”.
 - At the end of each mixing step, visually confirm that the solution is homogeneously mixed.

8. Procedure
The Fungitell STAT® product contains a Quick Visual Guide with illustrations and a summary of the features of the PKF08 instrument and BG Analytics® software.

The following procedures are already preset when using the PKF08 device and the BG Analytics® software: Device setting, evaluation of results and quality control. For more information, see the BG Analytics® Software User Manual or contact the manufacturer.

8.1 Instrument setting and test programming
8.1.1 When using PKF08 with BG Analytics® Software:
Turn on the device and follow the instructions of the BG Analytics® software. For detailed information, see the BG Analytics® manual.

- 8.1.2 When using another instrument and software, the following conditions should be met:
 - The instrument should be able to achieve and hold a temperature of 37°C±1°C.
 - The instrument and software must be able to read optical density over time (kinetic mode) at two wavelengths. Specifically, these wavelengths should be set to 405 nm and 495 nm.
 - Set the kinetic mode to a read length of 40 minutes (2400 seconds). Set the kinetic read interval to the minimum allowed by the software/instrument.
 - The measurement should be initiated immediately upon sample insertion.
 - Refer to the software manual to determine how to calculate a rate (slope) measurement from the data set. For the purposes of this test, this is generally achieved by executing a linear regression on the kinetic data over the time frame suggested. Set the linear regression calculation to execute over the range between 1900 and 2400 seconds using the “slice” function of the software.


- 8.2 Label tubes**
 - Label one empty tube for each patient serum sample to be tested.
 - Label one Fungitell STAT® Reagent tube for each patient serum sample to be tested.
 - Label one Fungitell STAT® Reagent tube for the Fungitell STAT® Standard.

- 8.3 Prepare patient serum sample**
 - Vortex patient serum samples for at least 20 seconds to ensure homogeneity.
Note: The freezing process can produce sample heterogeneity due to water abstraction to the growing ice crystal, thus excluding solutes.
 - Add the patient serum sample and Alkaline Pretreatment Solution in a ratio of 1:4 in the corresponding labeled empty tube. The recommended volumes are 50 µl of patient sample and 200 µl of Alkaline Pretreatment Solution.
Note: The Alkaline Pretreatment Solution converts triple-helix glucans into single-stranded glucans^{14,15} which are more reactive in the assay. Additionally, the alkaline pH serves to inactivate serum proteases and inhibitors that can interfere with the assay²⁴.
 - Vortex for 15 seconds and cover.

- 8.4 Prepare Fungitell STAT® Standard**
Note: Each product (Fungitell STAT® Standard and Fungitell STAT® Reagent pair) is tested and released independently. Thus, it is important to use the Lot# volumes of reconstitution and Alkaline Pretreatment Solution. These can be found on the Fungitell STAT® Standard package label, on the Fungitell STAT® product Certificate of Analysis, and available on the ACC website. Recommendation: Before starting the test, write down this information on the supplied Quick Visual Guide.
 - Reconstitute one vial of the Fungitell STAT® Standard with the Lot# specific volume of LAL Reagent Water and vortex for 15 seconds.
 - Add the Lot# specific volume of Alkaline Pretreatment Solution.
 - Vortex for 15 seconds and cover.

8.5 Pretreatment Incubation in tube reader
Incubate the patient serum sample tubes (from Step 8.3) and the Fungitell STAT® Standard vial (from Step 8.4) for 10 minutes at 37°C.

Note: When using the PKF08 instrument, on inserting a tube into a well, an indicator turns from red to green. Push the tube fully in until the indicator turns green.

 Caution, the tubes are fragile. In case of penetration of shards of glass and liquids into a measuring station of the PKF08, contact Associates of Cape Cod, Inc. Technical Service.

- 8.6 Prepare Fungitell STAT® Reagent tubes**
 - Reconstitute each of the Fungitell STAT® Reagent vials (labeled in Step 8.2 above) with 300 µl of LAL Reagent Water.
 - Vortex gently for no more than 5 seconds.
Note: The Fungitell STAT® Reagent contains a number of active proteins required for the assay and it is recommended to gently handle the solution. A maximum setting of 2000 RPM is recommended for any vortex device. Do not over mix.
 - At the end of the pre-incubation:
 - Transfer 75 µl of each patient serum sample solution into its corresponding Fungitell STAT® Reagent tube.
 - Transfer 75 µl of Fungitell STAT® Standard into its corresponding Fungitell STAT® Reagent tube.
 - Vortex all tubes for no more than 5 seconds and cover.

- 8.7 Start the run**
 - Insert the tubes into tube reader while confirming that each one is in the intended well.
 - Start the kinetic reading for a period of 40minutes, at 37°C.

9 Calculate the results
9.1 Measuring Principle
The results of the Fungitell STAT® test should be used as an aid in the diagnosis of an invasive fungal infection. The standard rates of the patient sample and Fungitell STAT® are derived from the calculation of the slope (rate) between 1900 and 2400 from the delta OD 405 - 495 nm results. The results of the Fungitell STAT® index are obtained from the division of the slope of the patient sample by the slope of the Fungitell STAT® Standard (see Figure 2).




Figure 2. Example of Fungitell STAT® kinetic curves and data analysis
The region highlighted in grey is the area of the slope determination (1900 to 2400 seconds), the solid line is an example Patient sample (PS) and the dashed line is the Fungitell STAT® Standard (STD). The slope of the sample (i.e., 0.00022 OD/s) divided by the slope of the 80 pg/mL Fungitell STAT® Standard (i.e., 0.00016 OD/s) leads to an Index of 1.4 for the sample.

- 9.2 When using PKF08 with BG Analytics® software:**
 - The review of the quality criteria is carried out automatically by the software. The result is displayed in the final report.
 - For valid test runs, the BG Analytics® software determines an index value for each sample, or assigns a clear negative or positive result to the sample.
 - If the software shows any indications of invalid parameters in the results evaluation, follow the instructions in the BG Analytics® Software Manual.

9.3 When using other software:
Verify that all Quality Criteria are met.

10. Quality Control
Listed below are the quality control criteria for the Fungitell STAT® Standard and patient serum sample results, including examples of expected kinetic curve shapes. These quality control criteria were validated during the studies presented in the Performance Features section.

- For all well numbers,** confirm Fungitell STAT® Standard or Sample # assignment

- For the Fungitell STAT® Standard result,**
 - the correlation coefficient (r) must be ≥ 0.980 and
 - the slope must be within the expected slope range of 0.00010 – 0.00024 OD/second.
If the Fungitell STAT® Standard result does not meet criteria #1 and #2, the run is invalid and all samples must be re-prepared and tested.

- For all patient sample results do the following:**
 - Determine if the result may be outside the measurement range of the test**
 - The result is likely out-of-range on the **positive** side if:
 - The Y intercept is positive and
 - The kinetic curve passes 0.4 OD before 1000 seconds.
 - The result is likely out-of-range on the **negative** side if:
 - The kinetic curve is positive after 500 seconds and
 - Has an OD ≥ 0.00 and < 0.07 at the end of the test.

*If the Sample result meets both criteria for either the positive or negative out-of-range, the general QC criteria below do not need to be completed, and the index value should **not** be calculated. All out of range results on the positive side should be reported as “Positive” and all out of range results on the negative side should be reported as “Negative”.*

- If the above criteria do not apply, verify the general QC:**
 - the kinetic curve must be positive after 500 seconds,
 - the kinetic curve must have an OD ≥ 0.00 at the end of the test,
 - the slope must be numerically positive,
 - the correlation coefficient (r) must be ≥ 0.980 and
 - the kinetic curve must have an upward increasing curve shape consistent with examples presented in **Figure 3**.

If the Sample result does not meet general QC criteria #1, 3-5, the sample result is invalid and the sample has to be tested again. Alternatively, a different method should be used. If the Sample result does not meet QC criteria #2, this suggests that the sample signal is low. In that case, the user should carefully review the provided curve in context and determine the validity of the results based on the laboratory internal quality system.

